

# HugeMed Series VL3S

**Display**  
Horizontal and vertical rotation angle:  $\geq 120^\circ$   
HD full view screen, 6 function buttons, USB port, HDMI port for big display.

**Camera System**  
Photo and Video recording function, Easy for clinical teaching and training

**Suction Button**  
Fast Connection; Smooth Suction

**Handle**  
Ergonomic streamline design  
Anti-slip and improve hand held comfort

**Working Chanel**  
Inner diameter: 1.2mm/2.6mm/3.0mm available

**Quick Handlings**  
1 lock design  
Fast installation and disassembling  
Stay firm and efficient

**Angulation Control**  
Flexible angulation control  
Precise location

**Scope**  
True 1:1 Torque Ratio for Easy operation

**Insertion part**  
Outer diameter: 2.5mm/3.9mm/5.2mm/5.8mm  
Fracture resistant, smooth, flexible and durable

**Distal tip part**  
Dual diflections  
New CMOS Technology  
High brightness LED light source  
160 up 130 down

## Reference

	Size	M58	M52	M39	M28
<b>Scope Body</b>	Illumination	LED			
	Joint Outer Diameter	5.8mm	5.2mm	3.9mm	2.8mm
	Instrument Chanel Inner Diameter	3.0mm	2.6mm	1.2mm	0mm
	Working Lenth	650mm			
	Anti-Fog	Anti-fog when it's on, no need for preheating.			
Operating system \ Software System upgrades are free for life					
<b>Display</b>	Size	3.5inch/8inch			
	Storage	32GB SD Card			
	Shoot	Support			
	Video Record	Support			
	White Balance	Support			
	Freeze/Backwards	Support			
	Video Synchronous Out put	Support			

Note: \* M28/M39 is not applicable for endoscopic diagnosis and treatment accessories.  
\* M28/M39/M52/M58 is not suitable for high-frequency electrotherapy surgery and laser hemostasis.  
\* Continuous working time in more than three hours, but it depends on ambient temperature and battery condition.

## Advantages

### Portable

Light weight, Easy to operate, Saving time for EMS, Convenient for transportation and storage

### Multifunction

Available for the intubation of difficult airway, Inspection of airway and double chamber, Alveolar lavage, Biopsy, Treatment of pneumothorax, Stent Implantation, Removal of tracheal foreign bodies, Other treatments for the bronchi, Such as Local radiotherapy and chemotherapy, Medication Injection

### Economic

One screen for four scopes, Four working channel sizes for option

### Durable

Electronic camera system with LED, New technology for the photo processing: CMOS, Higher quality pictures, Last longer than traditional scope

### Application

Applied to the departments such as: ICU, NICU, Respiration, Anesthesia, Emergency etc.

## Accessories available: Trolley

