



Single-use Choledochoscope

CH Series

Single-use



Cost-effective



No Cross-infection



Always ready to go



Safer, clearer, and faster Choledochoscopy

The excellent assistant during PTCD, PTCS, PTHC, LCBDE, and so on.



CH Series Imaging Part

Digital format with multiple output ports



Specification	
Model	HUV-01
Image output	1024 x 768 pixels
Screen ratio	4:3
Light source control	0-7 grades brightness control
Storage capacity	Maximum 16GB for dynamic or static imaging
Language	Chinese, English, Spanish, Portuguese, French
Power	100-240V~ 50/60HZ
Border type	7 kinds

Specification

Model	O.D.	I.D.	Deflection	Field of View	Working Length	Direction of View
CH-M52	5.0mm	2.6mm	210° Up/ Down	120°	350mm	0°

Ordering information

Product	Item code
Single-use Choledochoscope	
CH-M52	2019B-PA002
Image processor	
HUV-01	8002B-PA006

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.



Shenzhen HugeMed Medical Technical Development Co., Ltd.
415, 416-1, 516-1, Building 2, No.1, Mawu Road, Baoan Community, Yuanshan Street,
Longgang District, Shenzhen, Guangdong, 518115, China

+86-755-22275833 +86-755-22275833 info@hugemed.net

For more information
please visit our website:
www.hugemed.net